

Histone H3 Rabbit Polyclonal Antibody

Catalog No: #38103

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Histone H3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB,IHC,IF
Species Reactivity	Hu
Specificity	The Histone H3 Rabbit Polyclonal Antibody detects endogenous Histone H3 protein.
Target Name	Histone H3
Other Names	H3 histone antibody; HIST1H3A antibody; Histone cluster 1; H3a antibody
Accession No.	Swiss-Prot#:P68431
Uniprot	P68431
GeneID	8350;8351;8352;8353;8354;8355;8356;8357;8358;8968;
SDS-PAGE MW	15-17kd
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:2000

Images



Western blot analysis of Hela cells, using #38103 diluted at 1:2,000.

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Note: This product is for in vitro research use only