

C1orf57 antibody

Catalog No: #22295

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	C1orf57 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 127 and 190 of Human C1orf57
Target Name	C1orf57
Other Names	FLJ11383; MGC13186; RP4-659I19.2; RP4-678E16.2
Accession No.	Swiss-Prot:Q9BSD7 Gene ID:84284
Uniprot	Q9BSD7
GeneID	84284;
Concentration	0.8mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 21kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: MOLT4
12% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only