

RALB antibody

Catalog No: #22306

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	RALB antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 182 of Human RALB
Target Name	RALB
Accession No.	Swiss-Prot:P11234Gene ID:5899
Uniprot	P11234
GeneID	5899;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 23kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
 A: Hep G2
 12% SDS PAGE
 Primary antibody diluted at 1: 1000

Background

This gene encodes a GTP-binding protein that belongs to the small GTPase superfamily and Ras family of proteins. GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors. [provided by RefSeq]

Note: This product is for in vitro research use only