

RNASE11 antibody

Catalog No: #38185

Package Size: #38185-1 50ul #38185-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RNASE11 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	IHC
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous level of total RNASE11 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human RNASE11.
Target Name	RNASE11
Other Names	RNASE11; C14orf6;
Accession No.	Swiss-Prot#: Q8TAA1 NCBI Gene ID: 122651
Uniprot	Q8TAA1
GeneID	122651;
SDS-PAGE MW	22kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

IHC 1:50 - 1:100

Images



Immunohistochemistry of paraffin-embedded human prostate using RNASE11 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse lung using RNASE11 at dilution of 1:100 (40x lens).

Background

Belongs to the pancreatic ribonuclease family.

Note: This product is for in vitro research use only