

## NUMBL antibody

Catalog No: #38379

Package Size: #38379-1 50ul #38379-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	NUMBL antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total NUMBL protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human NUMBL.
Target Name	NUMBL
Other Names	CAG3A; CTG3a; NBL; NUMB-R; NUMBLIKE; NUMBR; TNRC23;
Accession No.	Swiss-Prot#: Q9Y6R0NCBI Gene ID: 9253
Uniprot	Q9Y6R0
GeneID	9253;
SDS-PAGE MW	65kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of Mouse kidney cells, using NUMBL .

## Background

---

Numb-like (Numbl or Nbl) binds Notch and is an intracellular modulator of Notch signaling pathway. Numbl contains a phosphotyrosine binding (PTB) domain that could bind the Y682 motif present in the cytoplasmic tail of APP. Deduced molecular weight of Numbl is 72 kDa.

---

Note: This product is for in vitro research use only