

## FAM50A antibody

Catalog No: #22337

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	FAM50A antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 149 and 337 of Human FAM50A
Target Name	FAM50A
Accession No.	Swiss-Prot:Q14320Gene ID:9130
Uniprot	Q14320
GeneID	9130;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 40kd

Western blotting: 1:500-1:3000

## Images



Sample (30 ug of whole cell lysate)

A: 293T

B: H1299

10% SDS PAGE

Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only