

GM130 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C07267Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GM130 Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 850-892 1002 derived from human GM130
Conjugates	Cy3
Target Name	GM130
Other Names	GM130; Golgin subfamily A member 2; 130 kDa cis-Golgi matrix protein; GM130 autoantigen; Golgin-95; GOLGA2
Accession No.	Swiss-Prot#:Q08379NCBI Gene ID:2801
Uniprot	Q08379
GeneID	2801;
Excitation Emission	512,550nm 570,615nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Golgi auto-antigen; probably involved in maintaining cis-Golgi structure.

Note: This product is for in vitro research use only