

TMSB4X antibody

Catalog No: #38656



Package Size: #38656-1 50ul #38656-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TMSB4X antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total TMSB4X protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human TMSB4X.
Target Name	TMSB4X
Other Names	FX; TB4X; PTMB4; TMSB4;
Accession No.	Swiss-Prot#: P62328NCBI Gene ID: 7114
Uniprot	P62328
GeneID	7114;
SDS-PAGE MW	5kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Immunofluorescence analysis of HeLa cells using TMSB4X .
Blue: DAPI for nuclear staining.

Background

This gene encodes an actin sequestering protein which plays a role in regulation of actin polymerization. The protein is also involved in cell proliferation, migration, and differentiation. This gene escapes X inactivation and has a homolog on chromosome Y.

Note: This product is for in vitro research use only