

BXDC1 antibody

Catalog No: #22367

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	BXDC1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 187 of Human BXDC1
Target Name	BXDC1
Accession No.	Swiss-Prot:Q9H7B2Gene ID:84154
Uniprot	Q9H7B2
GeneID	84154;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 36kd

Western blotting: 1:500-1:3000

Images



Sample (30ug whole cell lysate)
A: H1299
10% SDS PAGE
Primary antibody diluted at 1: 3000

Note: This product is for in vitro research use only