

# ETV7 antibody

Catalog No: #38781



Package Size: #38781-1 50ul #38781-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ETV7 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total ETV7 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ETV7.
Target Name	ETV7
Other Names	TEL2; TELB; TEL-2;
Accession No.	Swiss-Prot#: Q9Y603NCBI Gene ID: 51513
Uniprot	Q9Y603
GeneID	51513;
SDS-PAGE MW	38kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ETV7 at 1:1000 dilution.

## Background

---

The protein encoded by this gene belongs to the ETS family of transcription factors, which is a large group of evolutionarily conserved transcriptional regulators that play an important role in a variety of cellular processes throughout development and differentiation, and are involved in oncogenesis as well. This protein is predominantly expressed in hematopoietic tissues. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene (PMID:11108721).

---

Note: This product is for in vitro research use only