

CIAPIN1 antibody

Catalog No: #38835

Package Size: #38835-1 50ul #38835-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CIAPIN1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total CIAPIN1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human CIAPIN1 (NP_064709.2).
Target Name	CIAPIN1
Other Names	CIAPIN1;Anamorsin;DRE2;PRO0915;anamorsin
Accession No.	Uniprot:Q6FI81GeneID:57019
Uniprot	Q6FI81
GeneID	57019
SDS-PAGE MW	39kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using CIAPIN1 antibody.



Western blot analysis of extracts of various cell lines, using CIAPIN1 antibody.

Background

CIAPIN1 is a cytokine-induced inhibitor of apoptosis with no relation to apoptosis regulatory molecules of the BCL2 (MIM 151430) or CASP (see MIM 147678) families. Expression of CIAPIN1 is dependent on growth factor stimulation (Shibayama et al., 2004 [PubMed 14970183]).[supplied by OMIM, Mar 2008]

Note: This product is for in vitro research use only