

DTYMK antibody

Catalog No: #38855



Package Size: #38855-1 50ul #38855-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DTYMK antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total DTYMK protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human DTYMK.
Target Name	DTYMK
Other Names	CDC8; TMPK; TYMK; PP3731;
Accession No.	Swiss-Prot#: P23919NCBI Gene ID: 1841
Uniprot	P23919
GeneID	1841;
SDS-PAGE MW	23kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using DTYMK at 1:1000 dilution.



Immunofluorescence analysis of A549 cells using DTYMK .
Blue: DAPI for nuclear staining.

Background

Catalyzes the conversion of dTMP to dTDP.

Note: This product is for in vitro research use only