

MLN antibody

Catalog No: #38868

Package Size: #38868-1 50ul #38868-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MLN antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total MLN protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human MLN (NP_002409.1).
Target Name	MLN
Other Names	MLN;motilin
Accession No.	Uniprot:P12872GeneID:4295
Uniprot	P12872
GeneID	4295
SDS-PAGE MW	13kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded human stomach using MLN antibody.



Western blot analysis of extracts of various cell lines, using MLN Antibody.



Immunohistochemistry of paraffin-embedded mouse brain using MLN antibody.



Immunohistochemistry of paraffin-embedded mouse heart using MLN antibody.

Background

This gene encodes a small peptide hormone that is secreted by cells of the small intestine to regulate gastrointestinal contractions and motility. Proteolytic processing of the secreted protein produces the mature peptide and a byproduct referred to as motilin-associated peptide (MAP). Three transcript variants encoding different preproprotein isoforms but the same mature peptide have been found for this gene.

Note: This product is for in vitro research use only