

CLEC3B antibody

Catalog No: #38897

Package Size: #38897-1 50ul #38897-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLEC3B antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total CLEC3B protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human CLEC3B/Tetranectin (NP_003269.2).
Target Name	CLEC3B
Other Names	CLEC3B;TN;TNA
Accession No.	Uniprot:P05452GeneID:7123
Uniprot	P05452
GeneID	7123
SDS-PAGE MW	Refer to figures
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

IF□1:10 - 1:100

Images



Immunofluorescence analysis of HeLa cells using CLEC3B/Tetranectin antibody.

Background

Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.

Note: This product is for in vitro research use only