

DHRS2 antibody

Catalog No: #38924



Package Size: #38924-1 50ul #38924-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	DHRS2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total DHRS2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human DHRS2 (NP_878912.1).
Target Name	DHRS2
Other Names	DHRS2;HEP27;SDR25C1
Accession No.	Uniprot:Q13268GeneID:10202
Uniprot	Q13268
GeneID	10202
SDS-PAGE MW	27KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:10 - 1:100

Images



Western blot analysis of extracts of various cell lines, using DHRS2 antibody.



Immunofluorescence analysis of HeLa cells using DHRS2 antibody.

## Background

This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members of this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics. Alternative promoter use and alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only