

CD34 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C07644Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD34 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Description	KLH conjugated synthetic peptide aa 210-260 385 derived from human CD34
Conjugates	Cy5
Target Name	CD34
Other Names	Hematopoietic progenitor cell antigen CD34; CD34
Accession No.	Swiss-Prot#:P28906NCBI Gene ID:947
Uniprot	P28906
GeneID	947;
Excitation Emission	625,650nm 670nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Possible adhesion molecule with a role in early hematopoiesis by mediating the attachment of stem cells to the bone marrow extracellular matrix or directly to stromal cells. Could act as a scaffold for the attachment of lineage specific glycans, allowing stem cells to bind to lectins expressed by stromal cells or other marrow components. Presents carbohydrate ligands to selectins.

Note: This product is for in vitro research use only