

# FGF14 antibody

Catalog No: #39028



Package Size: #39028-1 50ul #39028-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FGF14 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total FGF14 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human FGF14 (NP_787125.1).
Target Name	FGF14
Other Names	FGF14;FGF-14;FHF-4;FHF4;SCA27
Accession No.	Uniprot:Q92915GeneID:2259
Uniprot	Q92915
GeneID	2259
SDS-PAGE MW	28kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using FGF14 antibody.

## Background

---

The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. A mutation in this gene is associated with autosomal dominant cerebral ataxia. Alternatively spliced transcript variants have been found for this gene.

---

Note: This product is for in vitro research use only