

FRZB antibody

Catalog No: #39030



Package Size: #39030-1 50ul #39030-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FRZB antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total FRZB protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human FRZB (NP_001454.2).
Target Name	FRZB
Other Names	FRZB;FRE;FRITZ;FRP-3;FRZB-1;FRZB-PEN;FRZB1;FZRB;OS1;SFRP3;SRFP3;hFIZ
Accession No.	Uniprot:Q92765GeneID:2487
Uniprot	Q92765
GeneID	2487
SDS-PAGE MW	36kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse skeletal muscle, using FRZB antibody.

Background

The protein encoded by this gene is a secreted protein that is involved in the regulation of bone development. Defects in this gene are a cause of female-specific osteoarthritis (OA) susceptibility.

Note: This product is for in vitro research use only