

NUDC antibody

Catalog No: #39094



Package Size: #39094-1 50ul #39094-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NUDC antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total NUDC protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human NUDC (NP_006591.1).
Target Name	NUDC
Other Names	NUDC;HNUDC;MNUDC;NPD011
Accession No.	Uniprot:Q9Y266GeneID:10726
Uniprot	Q9Y266
GeneID	10726
SDS-PAGE MW	45kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NUDC antibody.

Background

This gene encodes a nuclear distribution protein that plays an essential role in mitosis and cytokinesis. The encoded protein is involved in spindle formation during mitosis and in microtubule organization during cytokinesis. Pseudogenes of this gene are found on chromosome 2.

Note: This product is for in vitro research use only