

HPS6 antibody

Catalog No: #22400

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	HPS6 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 405 and 456 (Q86YV9) of HPS6
Target Name	HPS6
Accession No.	Swiss-Prot:Q86YV9Gene ID:79803
Uniprot	Q86YV9
GeneID	79803;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 83kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: A431
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Background

This intronless gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. This protein interacts with Hermansky-Pudlak syndrome 5 protein. Mutations in this gene are associated with Hermansky-Pudlak

syndrome type 6. [provided by RefSeq]

Note: This product is for in vitro research use only