

PKD/PKCm(Phospho-Ser910) Antibody

Catalog No: #11096



Package Size: #11096-1 50ul #11096-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PKD/PKCm(Phospho-Ser910) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of PKD/PKCm only when phosphorylated at serine 910.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of serine 910 (R-V-S(p)-I-L) derived from Human PKD/PKCm.
Target Name	PKD/PKCm
Modification	Phospho
Other Names	KPCD1; PKC-mu; PKCM; PKD; PRKCM
Accession No.	Swiss-Prot: Q15139NCBI Protein: NP_002733.2
Uniprot	Q15139
GeneID	5587;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 115kd


Western blotting: 1:500~1:1000

Immunofluorescence: 1:100~1:200

Images



Western blot analysis of extracts from EGF-treated 293 and PMA-treated HUVEC cells using PKD/PKCμ(Phospho-Ser910) Antibody #11096.



Immunofluorescence staining of methanol-fixed HeLa cells using PKD/PKCμ(Phospho-Ser910) Antibody #11096.

Background

Converts transient diacylglycerol. (DAG) signals into prolonged physiological effects, downstream of PKC. Involved in resistance to oxidative stress through activation of NF-kappa-B.

Matthews SA, et al. J Biol Chem 1999 Sep; 274(37): 26543-26549

Brandlin I, et al. J Biol Chem 2002 Feb; 277(8): 6490-6496

Storz P, et al. Mol Pharmacol 2004 Oct; 66(4): 870-879

Note: This product is for in vitro research use only