

HCE antibody

Catalog No: #22405

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	HCE antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 40 and 256 of Human RNGTT
Target Name	HCE
Accession No.	Swiss-Prot:O60942Gene ID:8732
Uniprot	O60942
GeneID	8732;
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 69kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: HeLa
B: Hep G2
C: MOLT4
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only