

# TCP-1 theta antibody

Catalog No: #22412



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TCP-1 theta antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 111 and 337 of Human CCT8
Target Name	TCP-1 theta
Accession No.	Swiss-Prot:P50990Gene ID:10694
Uniprot	P50990
GeneID	10694;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 60kd  
Western blotting: 1:500-1:3000  
Immunohistochemistry: 1:100-1:250

## Images



Sample (30 ug of whole cell lysate)  
A: A431  
B: H1299  
C: Hela  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded OVCAR3 xenograft, using TCP-1 theta antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only