

Plexin B1 Antibody

Catalog No: #39234



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com



Description

Product Name	Plexin B1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by protein A/G column
Applications	WB;IHC
Species Reactivity	Human; Mouse; Rat
Specificity	Plexin B1 antibody detects endogenous levels of Plexin B1 protein. Not tested with other proteins
Immunogen Type	Peptide-KLH
Immunogen Description	A peptide corresponding to the N terminal region of human PLXNB1 protein
Target Name	Plexin B1
Concentration	1.0mg/ml
Formulation	Reconstitute in 100uL of sterile H2O with 50% glyceroloO 0.01% thimerosal.
Storage	Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles

Application Details

WB 1:500;
IHC 1:100-200

Images



Note: This product is for in vitro research use only