

## NIR1 antibody

Catalog No: #22424

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	NIR1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 912 and 974 of Human PITPNM3
Target Name	NIR1
Accession No.	Swiss-Prot:Q9BZ71Gene ID:83394
Uniprot	Q9BZ71
GeneID	83394;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 107kd

Western blotting: 1:500-1:3000

## Images



Sample(30 ug of whole cell lysate)  
A: 293T  
B: Raji  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1500

## Background

PITPNM3 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the Drosophila retinal degeneration B (rdgB) protein (Ocaka et al., 2005 [PubMed 15627748]).[supplied by OMIM]

---

Note: This product is for in vitro research use only