

RPB8 antibody

Catalog No: #22443



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	RPB8 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 120 of RPB8
Target Name	RPB8
Accession No.	Swiss-Prot:P52434Gene ID:5437
Uniprot	P52434
GeneID	5437;
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 17kd  
Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)  
A: A431  
12% SDS PAGE  
Primary antibody diluted at 1: 1000

Background

This gene encodes one of the essential subunits of RNA polymerase II that is shared by the other two eukaryotic DNA-directed RNA polymerases, I and III. [provided by RefSeq]

---

Note: This product is for in vitro research use only