

ZKSCAN3 Antibody

Catalog No: #39589



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZKSCAN3 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human ZKSCAN3 protein 175-440aa (BC006118).
Target Name	ZKSCAN3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

WB

Images



Note: This product is for in vitro research use only