

CRHSP-24 antibody

Catalog No: #22458



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CRHSP-24 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 122 of human CARHSP1
Target Name	CRHSP-24
Accession No.	Swiss-Prot:Q9Y2V2Gene ID:23589
Uniprot	Q9Y2V2
GeneID	23589;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 16kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: Hep G2
15% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded A549 xenograft, using CRHSP-24 antibody at 1: 100 dilution.

Background

May play a role in acinar cell metabolism.

Note: This product is for in vitro research use only