

ZNF564 Antibody

Catalog No: #39690



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZNF564 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Affinity purified |
| Applications | ELISA WB IHC |
| Species Reactivity | Hu Rt Ms |
| Immunogen Type | Recombinant Protein |
| Immunogen Description | raised against a recombinant human ZNF564 protein 113-404aa (BC065023). |
| Target Name | ZNF564 |
| Concentration | 1.0mg/ml |
| Formulation | PBS, pH 7.4 with 0.02% Sodium Azide |
| Storage | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal intestine showing nuclear staining with anti-ZNF564 antibody at a dilution of 1/100.



Note: This product is for in vitro research use only