

GIPC3 antibody

Catalog No: #22474

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GIPC3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 1 and 14 of Human GIPC3
Target Name	GIPC3
Accession No.	Swiss-Prot:Q8TF64Gene ID:126326
Uniprot	Q8TF64
GeneID	126326;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 34kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: H1299
10% SDS PAGE
Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only