

PODN Antibody

Catalog No: #39852



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	PODN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal PODN (1) antibody was raised against a recombinated human PODN protein 69-236aa (BC030608).
Target Name	PODN
Concentration	1.0mg/ml
Formulation	Liquid in 1% BSA oO 0.02% NaN3oO sterile H2O with 50% glycerol.
Storage	Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.


Application Details

ELISA 1:20000-1:80000 WB 1:500-1:2000 IHC 1:100-1:200

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney staining with anti-PODN antibody at a dilution of 1/100



---

Note: This product is for in vitro research use only