

PPM1K antibody

Catalog No: #22478



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PPM1K antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 334 of Human PPM1K
Target Name	PPM1K
Other Names	DKFZp667B084
Accession No.	Swiss-Prot:Q8N3J5Gene ID:152926
Uniprot	Q8N3J5
GeneID	152926;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 41kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250

Images



Western Blot analysis of PPM1K expression in transfected 293T cell line by PPM1K polyclonal antibody.
A: Non-transfected lysate.
B: PPM1K transfected lysate.
12% SDS PAGE
Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded endometrial OVCA, using PPM1K antibody at 1: 100 dilution.

Note: This product is for in vitro research use only