

TFEC antibody

Catalog No: #22480

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TFEC antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 1 and 13 of TFEC
Target Name	TFEC
Accession No.	Swiss-Prot:O14948Gene ID:22797
Uniprot	O14948
GeneID	22797;
Concentration	0.4mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 36kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: 293T

10% SDS PAGE

TFEC antibody diluted at 1: 1000

Note: This product is for in vitro research use only