

MPZL3 Antibody

Catalog No: #39995



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | MPZL3 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Affinity purified |
| Applications | ELISA WB IHC |
| Species Reactivity | Hu |
| Immunogen Type | Recombinant Protein |
| Immunogen Description | Rabbit polyclonal MPZL3 (1) antibody was raised against a recombinant human MPZL3 protein 38-122aa (BC113586). |
| Target Name | MPZL3 |
| Concentration | 1.0mg/ml |
| Formulation | PBS, pH 7.4 with 0.02% sodium azide. |
| Storage | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |

Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000 IHC: 1:100-1:500

Images



Note: This product is for in vitro research use only