

GFPT1 Rabbit mAb

Catalog No: #40024



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	GFPT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GFPT1
Target Name	GFPT1
Concentration	1.0mg/ml
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.


Application Details

WB 1:1000 - 1:5000;  
IHC 1:100 - 1:200

Images



Western blot analysis of GFPT1 expression in MCF7 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

---

Note: This product is for in vitro research use only