

## UQCRFS1 antibody

Catalog No: #22495

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	UQCRFS1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 176 of Human UQCRC1
Target Name	UQCRFS1
Accession No.	Swiss-Prot:P31930Gene ID:7384
Uniprot	P31930
GeneID	7384;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 53kd

Western blotting: 1:500-1:3000

## Images



Sample(30 ug of whole cell lysate)  
 A: Hep G2  
 7.5% SDS PAGE  
 Primary antibody diluted at 1: 1000

## Background

This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This protein may mediate formation of the complex between cytochromes c and c1.

---

Note: This product is for in vitro research use only