

PDE1A antibody

Catalog No: #22496



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PDE1A antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 212 and 520 of PDE1A
Target Name	PDE1A
Accession No.	Swiss-Prot:P54750Gene ID:5136
Uniprot	P54750
GeneID	5136;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 61kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:50-1:500
Immunofluorescence: 1:100-1:200


Images




Sample (30 ug of whole cell lysate)
A: H1299
B: Hela
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded MDAMB468 xenograft, using PDE1A antibody at 1: 500 dilution.



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using PDE1A antibody at 1: 200 dilution.

Note: This product is for in vitro research use only