

## RASSF4 Antibody

Catalog No: #40071

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RASSF4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RASSF4 protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human Ras association (RaGDS/AF-6) domain family member 4
Target Name	RASSF4
Other Names	AD037
Accession No.	Swiss-Prot:Q9H2L5?Gene Accssion:BC032593
Uniprot	Q9H2L5
GeneID	83937;
SDS-PAGE MW	37KD
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:200-1:1000

## Images



Gel: 8%SDS-PAGE  
Lysate: 40ug Hepg2 cellPrimary antibody: 1/200 dilution  
Secondary antibody dilution: 1/8000  
Exposure time: 5 minutes

## Background

The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several

transcript variants; however, most of the variants have not been fully described.

---

Note: This product is for in vitro research use only