

TBRG1 Antibody

Catalog No: #40137

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TBRG1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TBRG1 protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	TBRG1
Other Names	NIAM; TB-5
Accession No.	Swiss-Prot:Q96LV8Gene Accssion:BC032312
Uniprot	Q96LV8
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40137 at dilution 1/35.

Background

Acts as a growth inhibitor. Can activate p53/TP53, causes G1 arrest and collaborates with CDKN2A to restrict proliferation, but does not require either protein to inhibit DNA synthesis. Redistributes CDKN2A into the nucleoplasm. Involved in maintaining chromosomal stability.WB analysis.

Note: This product is for in vitro research use only