

COG4 Antibody

Catalog No: #46966

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	COG4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total COG4 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human COG4
Target Name	COG4
Other Names	COD1; CDG2J
Accession No.	Swiss-Prot#:Q9H9E3NCBI Gene ID:25839Gene Accssion:BC013347
Uniprot	Q9H9E3
GeneID	25839;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46966(COG4 Antibody) at dilution 1/30. (Original magnification: ?00)

Background

The protein encoded by this gene is a component of an oligomeric protein complex involved in the structure and function of the Golgi apparatus. Defects in this gene may be a cause of congenital disorder of glycosylation type IIj. Two transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only