

DCAF13 Antibody

Catalog No: #47024

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DCAF13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DCAF13 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human DCAF13
Target Name	DCAF13
Other Names	GM83; WDSOF1; HSPC064
Accession No.	Swiss-Prot#:Q9NV06NCBI Gene ID:25879Gene Accssion:NP_056235
Uniprot	Q9NV06
GeneID	25879;
Calculated MW	51 kDa
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Western blotting:1:200-1000

Images



Gel: 8%SDS-PAGE

Lysate: 40 µg, Lane: Human fetal liver tissue lysate

Primary antibody:DCAF13 Antibody at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute

Background

DCAF13 localizes to the granular component of the nucleolus, a region associated with late ribosomal RNA-processing events and assembly of ribosomal particles.Possible role in ribosomal RNA processing (By similarity). May function as a substrate receptor for CUL4-DDB1 E3

ubiquitin-protein ligase complex.

Note: This product is for in vitro research use only