

DCAF5 Antibody

Catalog No: #47027

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DCAF5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DCAF5 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human DCAF5
Target Name	DCAF5
Other Names	BCRG2; BCRP2; WDR22; D14S1461E
Accession No.	Swiss-Prot#:Q96JK2 NCBI Gene ID:8816Gene Accssion:NP_003852
Uniprot	Q96JK2
GeneID	8816;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47027(DCAF5 Antibody) at dilution 1/20. (Original magnification: 200)

Background

May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.

Note: This product is for in vitro research use only