

DDX19A Antibody

Catalog No: #47036

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DDX19A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DDX19A protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DDX19A
Target Name	DDX19A
Other Names	DDX19L; DDX19-DDX19L
Accession No.	Swiss-Prot#:Q9NUU7NCBI Gene ID:55308Gene Accssion:BC005162
Uniprot	Q9NUU7
GeneID	55308;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47036(DDX19A Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47036(DDX19A Antibody) at dilution 1/20. (Original magnification: ?00)

Background

ATP-dependent RNA helicase involved in mRNA export from the nucleus. Rather than unwinding RNA duplexes, DDX19 functions as a remodeler of ribonucleoprotein particles, whereby proteins bound to nuclear mRNA are dissociated and replaced by cytoplasmic mRNA binding proteins (By similarity).

Note: This product is for in vitro research use only