

DLEC1 Antibody

Catalog No: #47059

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DLEC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DLEC1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human DLEC1
Target Name	DLEC1
Other Names	F56; DLC1; DLC-1; CFAP81
Accession No.	Swiss-Prot#:Q9Y238 NCBI Gene ID:9940Gene Accssion:NP_031361
Uniprot	Q9Y238
GeneID	9940;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47059(DLEC1 Antibody) at dilution 1/25. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 47059(DLEC1 Antibody) at dilution 1/25. (Original magnification: ?00)

Background

The cytogenetic location of this gene is 3p21.3, and it is located in a region that is commonly deleted in a variety of malignancies. Down-regulation of this gene has been observed in several human cancers including lung, esophageal, renal tumors, and head and neck squamous cell carcinoma. In some cases, reduced expression of this gene in tumor cells is a result of aberrant promoter methylation. Several alternatively spliced transcripts have been observed that contain disrupted coding regions and likely encode nonfunctional proteins.

Note: This product is for in vitro research use only