

DPT Antibody

Catalog No: #47077

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DPT Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DPT protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DPT
Target Name	DPT
Other Names	TRAMP
Accession No.	Swiss-Prot#:Q07507NCBI Gene ID:1805Gene Accssion:BC033736
Uniprot	Q07507
GeneID	1805;
Concentration	0.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47077(DPT Antibody) at dilution 1/20. (Original magnification: ?00)

Background

Dermatopontin is an extracellular matrix protein with possible functions in cell-matrix interactions and matrix assembly. The protein is found in various tissues and many of its tyrosine residues are sulphated. Dermatopontin is postulated to modify the behavior of TGF-beta through interaction with decorin.?

Note: This product is for in vitro research use only