

FAM198A Antibody

Catalog No: #47119

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FAM198A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FAM198A protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human FAM198A
Target Name	FAM198A
Other Names	C3orf41
Accession No.	Swiss-Prot#:Q9UFP1 NCBI Gene ID:729085Gene Accssion:NP_001123380
Uniprot	Q9UFP1
GeneID	729085;
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47119(FAM198A Antibody) at dilution 1/25. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47119(FAM198A Antibody) at dilution 1/25. (Original magnification: ?00)

Background

Note: LOC26172 (GeneID: 26172) was previously annotated as C3orf41, and a chimeric/readthrough mRNA that contained part of this gene, GeneID: 729085 and GeneID: 26172 was associated with the locus. The two genes are now represented by separate entries, LOC26172 (GeneID: 26172) and C3orf41 (GeneID: 729085).

Note: This product is for in vitro research use only