

## FAM91A1 Antibody

Catalog No: #47120

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FAM91A1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FAM91A1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human FAM91A1
Target Name	FAM91A1
Accession No.	Swiss-Prot#:Q658Y4 NCBI Gene ID:157769Gene Accssion:NP_659400
Uniprot	Q658Y4
GeneID	157769;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C

## Application Details

Immunofluorescence: 1: 20-100

## Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47120(FAM91A1 Antibody) at dilution 1/35. (Original magnification: ?00)

## Background

FAM91A1 (family with sequence similarity 91, member A1) is a 94 kDa phospho-protein that is not yet characterized.

---

Note: This product is for in vitro research use only