

## GPBP1L1 Antibody

Catalog No: #47129

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                  |
|-----------------------|------------------------------------------------------------------|
| Product Name          | GPBP1L1 Antibody                                                 |
| Host Species          | Rabbit                                                           |
| Clonality             | Polyclonal                                                       |
| Purification          | Antigen affinity purification                                    |
| Applications          | IHC                                                              |
| Species Reactivity    | Hu                                                               |
| Specificity           | The antibody detects endogenous levels of total GPBP1L1 protein. |
| Immunogen Type        | protein                                                          |
| Immunogen Description | Fusion protein of human GPBP1L1                                  |
| Target Name           | GPBP1L1                                                          |
| Other Names           | SP192                                                            |
| Accession No.         | Swiss-Prot#:Q9HC44NCBI Gene ID:60313Gene Accssion:BC015713       |
| Uniprot               | Q9HC44                                                           |
| GeneID                | 60313;                                                           |
| Concentration         | 0.7mg/ml                                                         |
| Formulation           | Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.  |
| Storage               | Store at -20C                                                    |


## Application Details

Immunofluorescence:1: 20-100

## Images



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47129(GPBP1L1 Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 47129(GPBP1L1 Antibody) at dilution 1/20. (Original magnification: ?00)

## Background

Possible transcription factor.

Note: This product is for in vitro research use only