

## MLLT6 Antibody

Catalog No: #47155

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MLLT6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MLLT6 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human MLLT6
Target Name	MLLT6
Other Names	AF17
Accession No.	Swiss-Prot#:P55198NCBI Gene ID:4302Gene Accssion:BC007237
Uniprot	P55198
GeneID	4302;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C

## Application Details

Immunofluorescence:1: 20-100

## Images



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47155(MLLT6 Antibody) at dilution 1/40. (Original magnification: ?00)

## Background

AF17 was discovered as part of a genetic rearrangement in the MLL (mixed lineage leukemia) gene associated with acute myeloid leukemia. AF17 is also known as MLLT6.

---

Note: This product is for in vitro research use only