

RADIL Antibody

Catalog No: #47185

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RADIL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RADIL protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human RADIL
Target Name	RADIL
Other Names	RASIP2
Accession No.	Swiss-Prot#:Q96JH8NCBI Gene ID:55698Gene Accssion:BC004919
Uniprot	Q96JH8
GeneID	55698;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 40-200

Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47185(RADIL Antibody) at dilution 1/40. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47185(RADIL Antibody) at dilution 1/40. (Original magnification: ?00)

Background

Downstream effector of Rap required for cell adhesion and migration of neural crest precursors during development.

Note: This product is for in vitro research use only